

Datasheet

V1.0
Product name: SARS-CoV-2 Spike Receptor Binding Domain (Spike-RBD) protein, C-terminal His tag

Catalog number	nCoV-p007	Quantity	100µg/200µg/500µg/1mg
Protein Information			
Origin	SARS-CoV-2	Buffer	PBS, pH=7.2
Expression Host	Expi293	Purity	≥85%
Theoretical MW	23kDa	Endotoxin	<10EU/mg
		Tag	C-His tag
Application	WB, ELISA, etc. Notice: please test in small quantity ahead of important experiment.		
Storage	2 - 8°C for up to 1 weeks. -20°C or even lower temperature for long term preservation. Avoid repeated freeze thaw cycles.		
Biosafety	This product is purified reconstructed protein with its original bioactivity. It has been filtered through 0.22µm filter membrane and does not contain any pathogenic bacteria.		
Important Notice	This product is for research use only. Please operate it under research lab protocol.		

Test data

Tris-Bis PAGE	
	