

Datasheet

V1.0

Product name: Human Angiotensin-Converting Enzyme 2 (ACE2) protein, C-terminal His tag

Catalog number	nCoV-P010	Quantity	100µg/200µg/500µg/1mg
Protein Information			
Origin	Human	Buffer	PBS, pH=7.2
Expression Host	Expi293	Purity	≥95%
Theoretical MW	77kDa	Endotoxin	<10EU/mg
Application	WB, ELISA, etc. Notice: please test in small quantity ahead of important experiment.		
Storage	2 - 8°C for up to 1 weeks. -20°C or even lower temperature for long term preservation. Avoid repeated freeze thaw cycles.		
Biosafety	This product is purified reconstructed protein with its original bioactivity. It has been filtered through 0.22µm filter membrane and does not contain any pathogenic bacteria.		
Important Notice	This product is for research use only. Please operate it under research lab protocol.		

Test data

